

IHE Integration Statement		Date	28-Dec-07
Vendor	Product Name	Version	
Ziosoft, Inc.	ZIOSTATION	V1.2	
This product implements all transactions required in the IHE Technical Framework (Revision 8.0) to support the IHE Integration Profiles Actors and Options listed below:			
Integration Profiles Implemented	Actors Implemented	Options Implemented	
Portable Data for Imaging	Portable Media Creator	None	
	Image Display	None	
	Print Composer	None	
Consistent Time	Time Client	None	
Internet address for vendor's IHE information: www.zio.co.jp/ihe			
Links to Standards Conformance Statements for the Implementation			
HL7	None		
DICOM	www.zio.co.jp/dicom		
Links to general information on IHE			
In North America: www.ihe.net	In Europe: www.ihe-europe.org	In Japan: www.ihe-j.org	