

<b>IHE Integration Statement</b>		<b>Date</b>	20-Sep-18
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	
Ziosoft, Inc.	Ziostation2 Plus	V2.9	
This product implements all transactions required in the IHE Technical Framework (Revision 14.0) to support the IHE Integration Profiles Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Portable Data for Imaging	Portable Media Creator	None	
	Image Display	None	
	Print Composer	None	
Consistent Time	Time Client	None	
<b><u>Internet address for vendor's IHE information:</u></b> <a href="http://www.zio.co.jp/ihe">www.zio.co.jp/ihe</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7</b>	None		
<b>DICOM</b>	<a href="http://www.zio.co.jp/products/dicom">www.zio.co.jp/products/dicom</a>		
<b>Links to general information on IHE</b>			
In North America: <a href="http://www.ihe.net">www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>	